

SUPAPORE - MPG



NEW PRODUCT RELEASE

SupaPore MPG are designed for general purpose use wherever a cost effective membrane filter is required. Constructed entirely from polypropylene it is an ideal cartridge for chemical and solvent resist applications in the Pharmaceutical and general Fine Chemical industries. Polypropylene membranes are also suitable for venting applications where a cost effective hydrophobic membrane is required.

Key Features

- High chemical resistance—all polypropylene construction
- Pleated membrane media for maximum efficiency
- High surface area for high flow and low clean pressure losses
- All sizes from 10 to 40" in modular construction
- Absolute ratings 0.1 and 0.2 μ m
- Cost effective 'versus' PTFE membrane
- Thermal bonded construction eliminates adhesives and ensures minimal extractables
- Identification on every cartridge for full traceability
- FDA compliant materials
- Available with all standard end cap configurations and seals

Technical Data

- Typical flow rates: 0.1 μ m 3 lt/min/69mBar
0.2 μ m 11 lt/min/69mBar
- Pre flushed with 17 megohm-cm water
- Hot water sanitisable @ 90°C
- Maximum pressure loss: Forward 3.4 bar @ 20°C
Reverse 2.7 bar @ 20°C

For further details please call the Sales Office on +44 (0)1276 670 600.

AMAZON FILTERS LTD.

ALBANY PARK ESTATE, FRIMLEY ROAD,
CAMBERLEY, SURREY, GU16 7PG
ENGLAND

Tel: +44 (0) 1276 670600
Fax: +44 (0) 1276 670101

Email: sales@amazonfilters.co.uk

Web: www.amazonfilters.co.uk



Cert No. FM39584